



# Gerbinol® CF

Tannin adsorption based on gelatines, silicates and isinglass, casein free


## Product description

Gerbinol® CF is a powdered, casein-free tannin adsorption agent. Gerbinol® CF evens out irregularities, rough edges and impurities. Gerbinol® CF alleviates and/or removes undesirable aromas and flavours. Gerbinol® CF is based on a variety of gelatines, silicates and isinglass.

Permitted according to EU Commission Regulation no. 934/2019. User must check compliance with national regulations. Laboratory tested for purity and quality.

## Dosage and use

Intended purpose	Dosage
Slight corrections	5 - 10 g/100 L
Tannin reduction	15 - 30 g/100 L
Correction of off-flavours	30 - 50 g/100 L



Gerbinol® CF can basically be used at any stage, i.e. including before conventional clarification processes. We recommend that treatment be carried out after clarification, however, as it is only then that flavour imperfections will be clearly discernable. Preliminary tests are recommended.

After the optimum dosage has been determined, Gerbinol® CF is stirred into ten times the volume of water (hand-warm if possible), mixing thoroughly so there are no lumps. Gerbinol® CF is then added to the main tank, stirring constantly. Stirring and dissolving of Gerbinol® CF should be performed with the care and diligence befitting a special product of this type.

## Storage

Store in a dry place. Protect from odours and moisture. Always close opened packs tightly.