



# e.Bois® Bags & Tubes

Pure, natural cotton bag for filling with e.bois® chips

## Product description

e.Bois® Bags and e.Bois® Tubes are pure cotton bags for filling with e.Bois® chips and then simply suspending in tanks, barrels and barriques. The aim is an easy contact between e.Bois® chips and wine and also very easy removal of e.Bois® chips from vessels.

## Notes on use

- The bags (for tanks) and tubes (for barrels) are made of washed and ironed pure cotton.
- They are easily wettable and completely reusable.
- They are very easy to use.
- After suspending they are completely immersed in the wines after just a relatively short time.
- It is easy to mix different qualities of e.Bois® in the ratios of your choice.
- e.Bois® Bags and e.Bois® Tubes are flexible, i.e. there is no restriction on barrel or bung diameter.
- The bags and tubes should be washed in hand-warm water again before first use.
- For easier emptying the tubes can be opened at both ends.
- The bags and tubes can be reused as many times as desired. We recommend washing them in hot water after each use. e.Bois® Bags and e.Bois® Tubes are boil-proof and so can be cleaned in the same way as normal laundry. This also precludes microbiological problems.
- The bags and tubes are suitable for use in washing machines. Detergent residues can be removed by washing several times in clean water.
- Bags and tubes used in red wine should not be used subsequently for white wine (risk of pigment residue transfer)
- The used chips can be reused as a decorative accessory, to cover soil in the garden or plant pots.
- The bags can also be used for filtering small quantities of sediment and as yeast bags.

## Dosage and use

- The bags can be filled with an arbitrary quantity of chips up to a maximum of 3 kg.
- The tubes can take up to 350 g chips; the easiest way to fill them is with a dry funnel.
- If greater quantities are required, several bags or tubes can be used in the same wine vessel.

## Storage

Away from odours. Wash and dry after use and if stored for extended periods.